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bre.flex 2nd Edition

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Instructions for use

bre.flex 2nd Edition

Before using the product, please read this product information carefully. Sold only to dental professionals

Why bre.flex 2nd Edition partial dentures?

No metal clasps – but gigiva-colored (pink) clasps that blend with the natural teeth.

- metal-free attachments (retention elements) protect the natural enamel
- better shade stability than comparable flexible base materials
- lasting plasticity and avoidance of embrittlement
- comfortable, resistant and durable
- entirely free from residual monomers and dimethacrylates. Perfectly suited for sensitive patients.

1. Indication range

bre.flex 2nd Edition is a thermoplastic material ans supplied as a granular material. The product can only be processed in the thermopress 400 injection molding system.

bre.flex 2nd Edition material is used to produce flexible partial dentures and unbreakable occlusal appliances. A range of four different colors (transparent, pink 2, pink 4 and pink-veined) is available to permit the fabrication of gingiva-colored mouldings.

2. Contraindication

Full dentures, crowns and bridges, attachments

2.1 Wax set-up

Use only wax for the set-up. Do not use acrylic or shellac base plates to avoid excessive thickness. bre.flex 2nd Edition material will not form a chemical bond with acrylic teeth. Therefore retentions need to be prepared in the acrylic teeth. Corresponding holes are drilled from the mesial to the distal direction using a Multidrill drill.

Fig. 1: 0.8 - 1.5 mm retention hole

 Fig. 2 : from mesial to distal

2.2 Wax model

Compared to conventional acrylic partial dentures, thinner are modelled for flexible partial dentures made of bre.flex 2nd Edition. The layer should just be thick enough to ensure sufficient stability especially in the area of critical junctions and the base. If properly designed and correctly processed, flexible partial dentures made of bre.flex 2nd Edition are unbreakable.

2.3 Recommended thickness of the wax model:

- Lingual thickness: 1.25 – 2.0 mm
- Palatal thickness: 0 – 1.5 mm
- Clasps: 1.0 – 1.5 mm (always uniform with margin and junction)
- Lingual margin: 1.25 – 2.0 mm
- Base: 0.8 – 1.75 mm
- Wax pattern should rest on the cingulum of anterior teeth.
- Wax pattern should include the posterior teeth sufficiently.

2.4 Attaching the injection channels (sprues)

Attach the sprues (v-shaped) using a 3.5 to 5 mm Biotec wax pattern (REF 430 0805 5; REF 430 0805 0) and connect them to form a main channel (sprue) (10 mm; resin sprue wax REF 430 0172 0).

Vaseline must be applied to the inside of the aluminium flask before the model is invested. Place the model in the stone as close to the injection hole of the flask as possible (stones of class III – IV). Stone quantity: approx. 350 g.

Fig. 03


Fig. 04

Fig. 05

2.5 Overbedding

Undercuts are carefully reduced using a cross-cut plaster bur. Apply a single coat of plaster against plaster insulating liquid to the plaster surface, close the flask and fill 350 – 400 g of bre.flex 2nd Edition into the flask. Stirring under vacuum and the use of a vibrating plate is recommended to avoid the formation of bubbles.

2.6 Boiling out

- Observe the setting time of the plaster in use (Thixo-Rock: 60 min)
- All four screws must be removed from the flask before boiling out!
- Heat at 90 °C for 7 – 8 min

- Open the flask carefully and boil out using clear, boiling water

All sharp edges must be smoothed off before the plaster surfaces are separated. The inner and outer sur-faces of the flasks sur must be contaminated with any wax.

2.7 Insulating

For insulating (separating) plaster against resin (bre.flex 2nd Edition) no alginate-based insulating liquids may be used. bredent offers two insulating liquids for the thermoplastic material bre.flex 2nd Edition:

- Apply a thin coat of Acrylic Sep (REF 520 0029 4) to the plaster surface which is still hot; resin bases with a silky matt luster are obtained.

- The varnish, transparent (REF 540 0100 6) is applied two times. Resin surfaces with a high luster are obtained.

3. Processing parameters

Program	Material	Temperature	Preheating time	Holding time	Speed	Force
5 Flex-2	bre.flex 2nd Edition	280 – 300° C	15 min	90 sec	6	165

- Preheat the thermopress 400 injection moulding system
- If the set temperature is reached, apply a thin coat of the thermopaste on the cartridge and insert into the heating chamber with the printing facing downward (2 – 5 mm deep into the chamber)
- Remove excess thermopaste.

Close the flask just 3 minutes before the end of the preheating time. Tighten the screws diagonally.

4. Deflasking

Deflasking should only be carried out using the devesting aid and the punch under the spindle press. This way careful handling and extended service life of the special flasks are ensured. All screws are removed prior to devesting. Press the material out of the flask halves under a spindle press using the devesting aid and the punch. Do not use a hammer! Remove the plaster from the vents in the upper part of the flask us-ing the flask hook and a hammer. If excess silicone grease is in the injection channel at the flask, it must be removed to avoid that these residues are injected into the cavity together with the material during the next injection process.

The injected dental restoration is devested with a large separating disc and plaster residues are removed with Dentaclean plaster solvent Speed (REF 520 0101 0).

4.1 Finishing and polishing

When cutting off the sprues, it must be ensured that the material will not become too hot. Work with high speed but do not exert any pressure. Use cross-cut bredent burs and diamond tools. bre.flex 2nd Edition can be perfectly processed with silicone burs by bredent (REF S263 QG 60). When polishing, do not exert pressure and use soft brushes, e.g. leather buff (REF 350 0099 0). Prepolishing with pumice powder and high luster polishing with Abraso-Star K50 (REF 520 0016 1) are performed.

5. Cleaning

When working with heated (ultrasonic) cleaning units it must be ensured that the temperature of the clean-ing liquid does not exceed 70° C. Do not use strong detergents! Do not clean with water steam (steam blaster).

Patient information

bre.flex 2nd Edition partial dentures are fabricated from a biocompatible nylon thermoplastic material. As a result, taking care of your new restoration is very easy. In order to maintain the perfect condition and natu-ral look of your restoration, the following recommendations should be observed:

Continue to brush your natural teeth and gums as directed by your dentist. Clean your restoration daily.

Simply rinse your restoration under warm running water for 30 – 60 seconds and soak it in bre.flex 2nd Edition cleaning bath for approx. 30 min.

Use only water and a soft toothbrush to clean your restoration. Do not use any whitening toothpaste or similar products since these may damage the surface of the material. When you are not wearing your res-toration, keep it in clean water. Please remember that although the denture base is unbreakable, the teeth are not! Please handle your restoration with care.

6. Safety recommendations and hazard warnings for dental technicians

6.1 Protective clothes: When processing bre.flex 2nd Edition in the thermopress injection moulding system, gloves must be worn. Risk of scalding!

6.2 Do not inhale vapors formed during preheating.

6.3 Caution: risk of formation of gas when overheating / or prolonged residence in the heating cylindr.

6.4 For safe handling of bre.flex 2nd Edition we refer to the material safety data sheet.

6.5 General safety precautions during processing:

In consideration of the maximally permissible residence time in the heating cylinder the maximum tempera-ture of the molten material should not be above 240° C for 5 min. If exposed to excess thermal and tempo-ral stress bre.flex will turn brown! Additionally, the pressure in the sealed cartridge may increase to such an extent that the lid is „blown off“ and persons may be harmed.

bre.flex 2nd Edition may not be mixed with other thermoplastic materials!

7. Storage and durability

Open package must be protected from moisture. If clean and dry storage is ensured, the granular material can be stored for an indefinite period. In a standard climate the absorption of water proceeds very slowly. Moist granular material is dried 4 hours at 80° C in the vacuum drying oven. Material must be used up with-in one year when package is opened.

8. Additional information

The information contained in these instructions for use is always updated according to the latest knowledge and experience. Therefore we recommend reading the instructions for use again before using a new pack-age. These non-binding user recommendations are based on our own experience. The user himself is re-sponsible for processing the product. Liability for incorrect results shall be excluded since bredent does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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GmbH & Co.KG · Weissenhorner Str. 2 · 89250 Senden · Germany